ORTOCENTER

ORTHOPEDIC PRODUCTS





























ORTOCENTER founded in 1985, one of the pioneers in the manufacture of ORTHOPEDIC and HOSPITAL PRODUCTS.

Company that has been highlighting the quality and technology applied in its products. Always researching and developing new techniques for evolution in the sector.

ORTOCENTER strictly maintains the sanitary surveillance standards required in the sector.

We have qualified professionals, equipment and machinery suitable for the production of quality products.

With these innovations new products have been developed expanding the range of options in the market.

We provide several rescue kits used in major emergency services nationwide.

Our quality has crossed borders and today we export worldwide.

We thank our customers, partners and employees for reliability, ensuring the quality and continuity of our products and company.

Director Celso R. Fieno

The Ortocenter caring for You

ORTOCENTER

ORTHOPEDIC PRODUCTS

Head and Neck

O3 Foam Cervical Collar
Cervical Collar Type Thomas
Cervical Collar with Mentonian Support
Emergency Cervical Collar
Philadelphia Cervical Collar
São Paulo Emergency Cervical Collar
Emergency Cervical Collar with Regulating

Trunk

- 04 Elastic Waistband with Duraluminium Fins
 Elastic Waistband with Flexible Fins
 Waistcoat Elastic Sirloin Putti Low
 Elastic Chest Band
 Thoracic Band in Tissue and Elastic
 Elastic Waist Band with Pendens
 Lumbar Band in Elastic with Crossed Straps
- 05 Putti High Elastic Vest
 Putti Vest in Denim
 Putti Vest in Denim and Spree Handle
 Vest Jewett
 Sprege for Adjustable Posture
 Denim Sprege
 Simple Sprege
 Abdominal Band 4 Panels

Upper limbs

- O6 Strap for Pregnant Woman
 Immobilizer in Eight for Clavicle
 Raimond Humeral Immobilizer
 Simple Sling Type Envelope
 Velpeau Type Immobilizer Nylon-Infant Sling
 Immobilizer Sling in Upholstered Brim
 Upholstered Nylon Immobilizer Sling
 Velpeau Type Immobilizer Sling on Canvas
- O7 Strap for Tennis Elbow
 Double Strap for Tennis Elbow
 Type Velpeau Immobilizing Sling on Canvas
 Sling Type Tote Bag on Canvas
 Sling Bag in Double Upholstered Denim
 Sling Type Tote Bag in Light Canvas-summer
 Tutor for Humeral Fracture
 Immobilizer Splint for Handle
- O8 Simple Elastic Wristband
 Long Immobilizing Splint
 Elastic Splint for Handle
 Ortocenter Elastic Splint
 Preventive Splint
 PVC Thumb Splint
 Hammer Finger Splint
 Splint Type Seesaw for Finger Extension

09 PVC Splint for Wrist and Thumb

Hands and Lower Limbs

09 PVC Splint for Wrist and Thumb
Elastic Hinged Knee Brace
Elastic Hinged knee Brace with Cross Straps
Partial Immobilizer Knee Orthosis
Immobilizing Orthesis Functional Position
PVC Splint for Hand and Finger Grip
PVC Splint for Handle
Immobilizing Splint Cuff and Thumb in Tissue

Lower Limbs

Appropriate to antique of the to the son of a

- 10 Ankle Stabilizer
 Ankle Immobilizer
 Elastic Compressive knee Brace
 Elastic Kneader for Kneecap
 Ankle Protector and Stabilizer
 Elastic Sub Patellar Strip with Cushion
 Sub Patellar Strip in Tube
 Plastic Tutor for Tibial Fracture
- 11 Anti Equine Support
 Damper for Calcaneal Spur
 Brace with Metatarsal Support
 Plantar Protector
 Sandal for Plaster
 Plantar Support
 Baruk Sandals
 Elastic Anklet in Eight
- 12 Boot Immobilizer Robocenter
 Short Immobilizer Boot
 Long Immobilizer Boot
 Elastic Strap for Traction of the Halux Valgo
 Gutter
 Insoles with Plantar Arches
 Insoles with Plantar Arches and Pilots
 or Metatarsal bar Retro Capital
 Center Insoles and Heels for Calcaneal Spur
- 13 Underwear for Inguinal Hernia Sling for Double Inguinal Hernia Sling for Inguinal Hernia Special Sling for Inguinal Hernia Sling for Umbilical Hernia Scrotal Suspenders

Ortoprene Line

Neoprene Coxal
Neoprene Elbow
Long Neoprene flat Knee Pad
Neoprene Hinged Knee Brace
with Cross Straps
Long Neoprene knee Pad with Holes
and Patellar U Support
Long Neoprene Anklet
Short Neoprene Anklet

ORTOCENTER

enjoyed of the freezeway is a margan

ORTHOPEDIC PRODUCTS

Line Neoprene

- 15 Neoprene wrist Protector with Thumb Guard
- 16 Lumbar Brace in Neoprene
 Neoprene Knee Pad with Tensioner
 Neoprene Glove with Thumb
 Neoprene Wrist Rest with Backing and
 PVC on Thumb
 Neoprene Non-thumb Munhequeira
 Simple Neoprene Thumb
 Neoprene Tennis Elbow Support
 Immobilizer Splice in Neoprene for Wrist

Child Line

17 Cervical Collar Schantz or Thomas
Cervical Collar with Mentonian Support
Structural Necklace
Brace for Umbilical Hernia
Frejka Diaper
Velpeau Type Immobilizer Nylon-Infant Sling
Pavlik Suspensory
Elastic Sling Handset for Lower
Limbs Rotation

Rescue and Emergency

- 18 Emergency Cervical Collar with Regulating
 Belt for Fastening on Plank Removal
 Dorsal Immobilizer Vest KED
 Dorsal Immobilizer Vest KED Child
 Emergency Plank in Naval Plywood
 Children's Polyethylene Emergency Board
 Long Plank for Immobilization in Polyethylene
 Emergency Sacolão
- 19 Pre-hospital Care Fellowship Hospital Care Backpack Medicine Backpack Trauma Bag

Silicone Line

20 Silicone Heel without Edge Silicone Heel Silicone Insole with Soft Bow Silicone Flat Insole Silicone Finger Separators

Several

21 Abdominal Brace for Obese and Post Surgical
Lumbar Brace with Suspensory Suspenders for
Work Safety
Thrat Necklace-Bath Protector
Brace for Neoprene Therapy
Brace for Pregnant Women with
Lumbar Reinforcement
Anti Snoring Track
Immobilizer Wrist Splint and Neoprene Thumb

Appen in in some of aminifficence in

22 Hip Abduction Pad
Egg Shell Mattress
Elbow Protector
Heel Protector
Bracelet and Anklet for Restrictions
Splice Type Multi Splint
Disposable Plus Restrictor

Rescue and Emergency

- 23 Triangular Bandage
 Immobilizer Head
 CIPA Type Immobilization KIT
 KIT Childbirth
 KIT Rescuer
 Aluminized Thermal Blanket
 Splice for Immobilization Aramada E.V.A.
 Splice for Immobilization Aramada Foam
- 24 Belts for Spider Type Surfboard Belts for Surfboard 3 Belts

השולם העול לבים מו בלבים מוחום בלי מוח ביות ביות

all second of the part of the second of the

profit find atraffigure trepoper to take for the gomes in the tree to forth some forthe

we whope to leave the terms to be by parter and from the lange for affects of these

along the sales wines the remain for the former of sales groups after more atternance of the sales

tale for I smarter atma from Johnson by the from the grave I tome Tome

Foam Cervical Collar



Cod: 005

Indications:

Traumatisms, torticollis, arthritis, arthroses and Somatizations.

Features:

Made of high density soft foam and Velcro closure.

Sizes: Extra Small, Small,

Medium, Large

Cervical Collar Schantz or Thomas



Cod: 010

Indications:

Traumatisms, torticollis, arthritis, arthroses and Somatizations.

Features

Made of high density soft foam and Velcro closure.

Sizes: Extra Small, Small,

Medium, Large

Cervical Collar with Mentonian Support



Cod: 015

Indications:

Traumatisms, torticollis, arthritis, including rheumatoid, arthroses and other cervical spine affections.

Features:

Made of Sturdy flexible plastic, upholstered on the edges, Mentonian support to reduce movement, height adjustment and velcro closure.

Sizes: Extra Small, Small, Medium, Large

Emergency Cervical Collar



Cod: 022

Indications:

Support for the cervical spine during removal of the injured traumatisms, torticollis, arthritis, arthroses and other cervical spine affections.

Features:

Made of Sturdy flexible plastic, has frontal orifice for tracheotomy procedure and velcro closure.

Sizes: Extra Small, Small, Medium, Large

Cervical Collar Philadelphia



Cod: 025

Indications:

Traumatisms, torticollis, arthritis, arthroses and other cervical spine affections.

Features:

Made of Polyformial foam by propositing control of bending, extension and rotation. Perforated for better ventilation and velcro closure.

Sizes: Baby, Extra Small, Small, Medium, Large

Emergency Cervical Collar



Cod: 035

Indications:

Support for the cervical spine during removal of the injured traumatisms.

Features:

Made of High density polyethylene, coated in E.V.A. and Mentonian support. Posui wide frontal opening for analysis of the choroid pulse, access to trachea with regulation in the circumference and velcro closure.

Sizes: Extra Small, Small, Medium, Large and Extra Large

Emergency Cervical Collar - Rescue



Cod: 037

Indications:

Support for the cervical spine during removal of the injured traumatisms.

Features:

Made of High density polyethylene, coated in E.V.A. and Mentonian support. With wide frontal opening for analysis of the choroid pulse, access to trachea and velcro closure.

Sizes: Extra Small, Small, Medium, Large and Extra Large

Emergency Cervical Collar - São Paulo



Cod: 041

Indications:

Support for the cervical spine during removal of the injured traumatisms. Traumatisms, torticollis, arthroses and other cervical spine affections.

Features:

Made of High density polyethylene, coated in E.V.A. 5mm for best confort. With wide frontal opening for analysis of the choroid pulse and velcro closure.

Sizes: Single with adjustable

Elastic Chest Band



Indications:

Contusions, rib fractures, steoporotic thoracic vertebrae.

Features:

Made of elastic and velcro closure.

Height: 15, 20 and 25 cm

Sizes: Small, Medium, Large

Cod: 055

Thoracic Band in Tissue and Elastic



Indications:

Contusions, rib fractures, steoporotic thoracic vertebrae.

Features:

Made of Linen Tissue, elastic and velcro closure.

Heights: 15, 20 and 25 cm

Sizes: Small, Medium, Large

Elastic Waist Band with Siddles and Lateral Reinforcement



Indications:

Abdominal wall sagging, post partum and post surgical.

Features:

Made of elastic With siddles, lateral reinforcement and velcro closure.

Heights: 20 cm

Sizes: Small, Medium, Large

and Extra large

Elastic Waist Band with Cross Straps



Cod: 063

Indications:

Abdominal wall sagging, post partum and post surgical.

Features:

Made of elastic With siddles, lateral reinforcement and velcro closure.

Heights: 25 cm

Sizes: Small, Medium, Large and

Extra large

Elastic Lumbar Range with Cross Straps



Indications:

Lumbar spine Stabilization, post partum and post surgical.

Features:

Made of elastic With siddles, lateral reinforcement and velcro closure.

Heights: 25 cm

Sizes: Small, Medium, Large

and Extra large

Elastic Waistcoat Strap with Duraluminium Fins



Indications:

Lumbar spine stabilization in traumatisms arthroses, scoliosis, fracture of other diseases of the region, stabilization of the trunk.

Features:

Made of elastic with siddles, lateral reinforcement, four duraluminium fins and velcro closure.

Tamanhos: Small, Medium, Large and Extra large.

Elastic Waistcoat Strap with Flexible Fins



Indications:

Lumbar spine stabilization in traumatisms arthroses, scoliosis, fracture of other diseases of the region, prevention of recurrence in sports practices.

Features:

Made of elastic with siddles, lateral reinforcement, auto adjustable fins with steel spring and velcro closure.

Sizes: Small, Medium, Large and Extra large.

Vest in Elastic Sirloin Sacrum Low Putti



Indications:

Stabilization of the sacral spine in case of contractures, post-trauma, arthroses, spondylolithesis, lordosis, postoperative and other diseases of the region.

Features:

Made of elastic with siddles, lateral reinforcement, duraluminium fins and velcro closure.

Sizes: Small, Medium, Large and Extra large.

Vest in Elastic Sirloin Sacrum High Putti



Cod: 074

Indications:

Stabilization of the region in case of herniated disc, post-trauma, arthroses, spondylolithesis, lordosis, postoperative and other diseases

Features:

Made of elastic with siddles, lateral reinforcement, duraluminium fins and velcro closure.

Small, Medium, Large and Extra large

Simple Backrest



Cod: 075

Indications:

Correction and posture of the shoulders and treatment of p ostural and juvenile thoracic kyphosis.

Features:

Made of sturdy elastic with foam auxiliary pads, Velcro closure

Small, Medium, Large and Extra large.

Backrest for Adjustable Posture



Cod: 078

Indications: Correction and posture of the shoulders and treatment of postural and juvenile thoracic kyphosis.

Features: Made of sturdy elastic with passant adjust, Velcro closure

Sizes: Small and Large

Jeans Backrest



Cod: 080

Indications:

Correction and posture of the shoulders and treatment of postural and juvenile thoracic kyphosis.

Features:

Made of sturdy elastic, jeans, Velcro closure.

Sizes:

Extra Small, Small, Medium, Large and Extra large.

Abdominal Band 4 Panels



Cod: 088

Post operative, abdominal liposuction, stomach reduction, hernias, post partum, bruising in the chest region and avoids sagging.

Features: Made of elastic panels, Velcro closure

Sizes: Small, Medium, Large and Extra large.

Heights 30 cm

Indications:

Stabilization of the region

Dorsum-Loin-sacrum in case of herniated disc, post-trauma, arthroses,

Putti Vest with Handle Type Sprege



Cod: 090

Indications:

Stabilization of the region Dorsum-Loin-sacrum in case of herniated disc, post-trauma, arthroses, spondylolithesis, lordosis, postoperative and other diseases of the region.

Features: Made of Jeans With side adjustments, steel fins, jean strap and Velcro closure.

Sizes: Small, Medium, Large and Extra large.

Putti Vest in Jeans



Cod: 100

Made of Jeans With side adjustments. steel fins and Velcro closure.

spondylolithesis, lordosis, postoperative and other diseases of the region.

Features:

Small, Medium, Large and Extra large.

Jewett Vest



Indications:

For patients with cyphotic posture and fracture in the thoracic spine.

Features:

Made of PVC and coated in soft foam and Velcro closure.

Sizes:

Small, Medium, Large and Extra large.

Cod: 105

Immobilizer in Eight for Clavicle



Código: 110

Indications:

Immobilization and alignment of clavicle fractures

Made with adjustable foam cotton strips and metal walkers with auxiliary pads.

Extra Small, Small, Medium and Large

Strap for Pregnant Woman



Código: 115

Indications:

Abdomen support, relief of low back pain, relief of abdominal muscles and spine. Posture Improvement.

Features:

Made with Front support in solid tissue with elastic rods.

Sizes:

Single Size

Simple Sling Type Envelope



Código: 200

Indications:

Upper Limb Support.

Features:

Made with Tissue.

Small, Medium and Large

Velpeau Type Immobilizer in Upholstered Nylon (Infant)



Código: 201

Indications:

Immobilization of the upper limb, contusions dislocation, fractures of the humerus, elbow. bursitis and sequelae of the A.V.C.

Features:

Made of nylon with foam, plastic buckle with one-handed manipulation, Velcro closure, bilateral.

Single Size - Child

Raymond Humerus Immobilizer



Código: 205

Indications:

Immobilize fractures of the humeral neck in the final phase of consolidation, dislocation and contusions of the shoulder, elbow and other affections of the upper limb. Stabilizes arm, forearm and shoulder

Features:

Upholstered tissue, velcro closure and elastic back

Sizes:

Small, Medium and Large

Velpeau Type Immobilizer in **Upholstered Tissue**



Código: 210

Indications:

Immobilization of the upper limb, contusions dislocation, fractures of the humerus, elbow, Dislocations, bursitis and sequelae of the A.V.C.

Features:

Made of upholstered tissue double, plastic buckle with one-handed manipulation, Velcro closure, bilateral.

Sizes:

Small, Medium and Large

Velpeau Type Immobilizer in **Upholstered Nylon**



Indications:

Immobilization of the upper limb, contusions dislocation, fractures of the humerus, elbow, Dislocations, bursitis and sequelae of the A.V.C.

Features:

Made of nylon, plastic buckle with one-handed manipulation, Velcro closure, bilateral.

Sizes:

Small, Medium and Large

Velpeau Type Immobilizer in Canvas



Código: 220

Indications:

Immobilization of the upper limb, contusions dislocation, fractures of the humerus, elbow, Dislocations, bursitis and sequelae of the A.V.C.

Features:

Made of Canvas, plastic buckle with one-handed manipulation, Velcro closure, bilateral.

Sizes:

Small, Medium and Large

Velpeau Type Immobilizer in Perforated Screen



Cod: 222

Indications:

Immobilization of the upper limb, contusions dislocation, fractures of the humerus, elbow, Dislocations, bursitis and sequelae of the A.V.C.

Features:

Made from perforated screen, plastic buckle with one-handed manipulation, Velcro closure, bilateral

Sizes: Small, Medium, Large and Extra Large

Sling Type Tote Bag in Perforated Screen



Cod: 225

Indications:

Upper limb support.

Features:

Made of Perforated screen, plastic buckle and velcro closure. Bilateral.

Sizes:

Small, Medium, Large and Extra Large

Sling Type Tote Bag in Double Upholstered Tissue



Cod: 226

Indications:

Upper limb support.

Features:

Made of Tissue, plastic buckle and velcro closure. Bilateral.

Sizes:

Small, Medium, Large and Extra Large

Velpeau Type Immobilizer in Canvas Summer Light Type



Cod: 230

Indications:

Upper limb support.

Features:

Made of Canvas, plastic buckle and velcro closure.
Bilateral.

Sizes

Small, Medium, Large and Extra Large

Strap for Tennis Elbow



Cod: 235

Indications:

Epicondylites and other affections of the elbow, prevention of recurrence of sports practice.

Features:

Upholstered Velcro strap with microfoam. Bilateral.

Sizes:

Small and Large

Double Strap for Tennis Elbow



Cod: 238

Indications:

Epicondylites and other affections of the elbow, prevention of recurrence of sports practice.

Features:

Strap with support on the handle and steel fins. Bilateral.

Sizes:

Small and Large

Tutor for Humeral Fractures



Cod: 250

Indications:

Functional treatment of humeral fractures by the Sarmiento method and assists in the stabilization of the humerus.

Features:

Made of semi rigid plastic molded in two pieces, mesh liner and Velcro closure.

Sizes:

Medium and Large. (Right and left)

Immobilizer Splint for Handle



Indications:

Tenosynovitis for wrist, arthroses, ligamentous lesions, corrections of the vicious positions in rheumatoid arthritis, region affections and prevention of recurrence in the return of activities.

Features:

Made of Tissue, Velcro closure and internal reinforcement in rigid plastic.

Sizes:

Small, Medium, Large and Extra large (Right and left)

80

Long Immobilizing Splint



Cod: 308

Indications:

Traumatic injuries, tenosynovitis for wrist, arthroses, ligamentous lesions, region affections and prevention of recurrence in the return of activities.

Features:

Made of Tissue, Velcro closure and internal reinforcement in rigid plastic.

Sizes:

Small, Medium, Large and Extra large (Right and left)

Elastic Splint for Handle



Cod: 310

Indications:

Traumatic injuries, tenosynovitis for wrist, region affections and prevention of recurrence in the return of activities.

Features:

Made of Elastic and Tissue, Velcro closure and internal reinforcement in rigid plastic.

Sizes:

Small, Medium, Large and Extra large. (Right and left)

Ortocenter Elastic Splint



Cod: 312

Indications:

Traumatic injuries, tenosynovitis, synovial cysts, region affections and prevention of recurrence in the return of activities.

Features:

Made of Elastic with support in Curvin, Velcro closure and internal reinforcement in rigid plastic.

Sizes:

Small, Medium, Large and Extra large (Right and left)

Preventive Splint



Cod: 315

Indications:

Prevention of Tendoinites during the work L.E.R.

Features:

Made of Tissue, Velcro closure and internal reinforcement in rigid plastic.

Sizes:

Small, Medium, Large and Extra large. (Right and left)

Immobilizer Cuff and Thumb Splint



Cod: 307

Indications:

Tenosynovitis for wrist, arthroses, ligamentous lesions, rheumatoid arthritis, region affections.

Features:

Made of Tissue, Velcro closure and internal reinforcement in rigid plastic.

Sizes:

Small, Medium, Large and Extra large. (Right and left)

Simple Elastic Wristband



Cod: 318

Indications:

Light affections of the wrist, in the sports practice prevention of recurrence.

Features:

Made of Elastic, Velcro closure.

Sizes:

Small, Medium, Large

Splint Type Seesaw for Finger Extension



Cod: 330

Indications:

Sprains, contusions, small fractures, joint reeducation in rheumatoid arthritis and in the arthroses and affections of the fingers.

Features:

Made of Stainless steel wire and micro porous foam pads.

Sizes:

Extra Small, Small, Medium, Large and Extra large

Hammer Finger Splint



Cod: 333

Indications:

Rupture of the extensor tendon (hammer finger) joint reduction of rheumatoid arthritis deformities.

Features:

Made of Stainless steel wire, Leather Support and E.V.A.

Sizes:

Extra Small, Small, Medium, Large and Extra large

PVC Splint for Thumb



Cod: 320

Indications:

Contusions, small fractures, joint reeducation in rheumatoid arthritis and in the arthroses and affections of the Thumb.

Features:

Molded PVC and Velcro tie Rods.

Single Size (Right and left)

PVC Splint for Handle and Thumb



Cod: 335

Indications:

Contusions, fractures, Tenosynovitis of Quervain, sprains, joint reeducation in rheumatoid arthritis and in the arthroses and affections of the Thumb.

Features:

Molded PVC and Velcro tie Rods.

Single Size (Right and left)

PVC Splint for Handle Hand and Thumb



Cod: 340

Indications:

Contusions, fractures, Peripheral nerve palsy Sequelae of A.V.C., joint reeducation in rheumatoid arthritis and in the arthroses and other affections. Finger and wrist immobilization.

Features:

Molded PVC and Velcro tie Rods.

Single Size (Right and left)

PVC Splint for Handle



Features:

Contusions, fractures on handle.

Indications:

Molded PVC and Velcro tie Rods.

Sizes:

Small, Medium and Large (Right and left)

Fixed functional position.

Partial Immobilizer Knee Orthosis



Traumatic knee disorders, prevention of recurrence in the return of activities. Reeducation of gait in sequelae of STROKE, arthroses, post-operative arthritis of cerebral palsy.

Features:

Made of Elastic, Tissue and foam, Velcro closure and internal reinforcement in Duralumin.

Small, Medium, Large

Immobilizing Orthosis Functional Position



Cod: 430

Indications:

Traumatic knee disorders, Gypsum replacement, prevention of recurrence in the return of activities, reeducation of gait in sequelae of STROKE, arthroses, post-operative arthritis of cerebral palsy.

Made of Elastic, Tissue and foam, Velcro closure and internal reinforcement in Duralumin.

Sizes:

50, 60, 70cm (height) - Narrow and wide

Elastic Hinged Knee Brace



Indications:

Traumatic knee disorders, Chronic instability, arthroses, arthritis.

Features:

Made of Elastic, side Duralumin fin and medial, Velcro closure.

Sizes:

Small, Medium, Large and Extra Large

Elastic Hinged Knee Brace with **Cross Straps**



Indications:

Traumatic knee disorders, Chronic instability, Recurrence Prevention of sports practice.

Features:

Made of Elastic, side Duralumin fin and medial, Crossstraps to enhance knee rotational stability, Velcro closure.

Sizes:

Small, Medium, Large and Extra Large

0

Elastic Compressive knee Brace



Indications:

Post-operative, post-arthroscopy, compression of joint effusion.

Features:

Made of Elastic, Hole for the kneecap and foam liner, Velcro closure.

Sizes:

Small, Medium, Large and Extra Large

Elastic Kneader for Kneecap



Cod: 455

Indications:

Chronic instability post-operative, post-arthroscopy.

Features:

Made of Elastic straps, Hole for the neecap and foam liner, Velcro closure.

Sizes:

Small, Medium, Large and Extra Large

Elastic Sub Patellar Strip with Cushion



Indications:

Tendinitis of the patellar tendon, stress relief in the patellar tendon.

Features:

Made of Elastic straps and micro foam, Velcro closure.

Sizes:

Small, Medium, Large

Sub Patellar Strip in Tube



Cod: 465

Indications:

Tendinitis of the patellar tendon, stress relief in the patellar tendon.

Features:

Made of Elastic straps with anterior cushion and micro foam, Velcro closure.

Sizes.

Small, Medium, Large and Extra large

Tutor for Tíbia Fractures



Indications: Functional tro

Functional treatment of Tíbia fractures by the Sarmiento method and Protection of the region after placement of prostheses or implants.

eatures:

Made of semi rigid plastic molded in two pieces, Articulated ankle support and Velcro closure.

Sizes:

Small and Large. (Right and left)

Ankle Stabilizer



Cod: 500

Indications:

Avoid sprains by improving the stability of the ankle joint in sports practice. Prevention of deformities of rheumatoid arthritis, assisting in the gait of the patient with flexor dorsum.

Features:

Made of Cotton canvas plastic side fins, elastic straps and Velcro closure.

Sizes:

Small, Medium, Large and Extra Large (Right and left)

Protect and Ankle Stabilizer



Avoid sprains by improving the stability of the ankle joint in sports practice.

Features:

Indications:

Made of Cotton canvas plastic side fins, elastic straps and Lock with Shoelace.

Sizes:

Small, Medium, Large and Extra Large. (Right and left)

Ankle Immobilizer



Cod: 510

Indications:

Mild sprains, replaces gypsum, prevention and treatment of rheumatoid arthritis deformities.

Features:

Made of leather, duraluminum side fins and Velcro closure.

Sizes:

Small, Medium, Large and Extra Large. (Right and left)

Elastic Anklet in Eight



Cod: 515

Indications:

Avoid sprains by improving the stability of the ankle joint in sports practice. Compression of the ankle joint, treatment of small sprains, tendinitis, bursitis.

Features:

Made of Elastic straps and Velcro closure.

Sizes: Single Size

Sandal for Plaster



Cod: 520

Indications:

Support and protection of plasted boots when allowed to march.

Features:

EVA sole, resistant synthetic tissue. Velcro closure.

Sizes:

Extra Small, Small, Medium, Large

Anti Equine Support



Cod: 526

Indications:

Reeducation in locomotion in AVC sequelae, Nerve paralysis.

Features:

Elastic, With Latex backing on the sole and velcro closure

Sizes:

Small, Medium, Large

Brace with Metatarsal Support



Cod: 530

Indications:

Anterior arch drop, hyperkeratosis in the metatarsal Vhead region and rheumatoid arthritis deformity.

Features:

Elastic, Latex Foam or silicone and velcro closure.

Sizes:

Small, Medium, Large

Plantar Support



Cod: 535

Indications:

Transverse arch Fall, calosities and ulcerations in the lumbar region of the metatarsal head, deformity in rheumatoid arthritis.


Features:

Latex Foam Coated on Jersey.

Sizes:

Small, Large

Protect Support



Cod: 540

Indications:

Transverse arch Fall, calosities and ulcerations in the lumbar region of the metatarsal head, deformity in rheumatoid arthritis.

Features:

Latex Foam Coated on Jersey.

Sizes:

Small, Large

Damper for Calcaneal Spur



Cod: 1545

Indications:

Plantar facite, heel bursitis, Calcaneal Spur.

Features:

Molded Micro Porous Foam.

Sizes:

Small, Medium, Large and Extra Large

Baruk Sandals



Cod: 550

Indications:

Post operative of Bunion.

Features:

Rubber, leather, and velcro closure.

Sizes

Extra Small, Small, Medium, Large

Insoles with Plantar Arches



Indications:

Pain relief and correction of foot deformities.

Features:

Micro Foam.

Sizes:

Numbers from 3 1/2 to 12 1/2

Insoles with Plantar Arches or Metatarsal Bar



Cod: 1810

Indications:

Pain relief and correction of foot deformities, metatarsalgia.


Features:

Micro Foam.

Sizes.

Numbers from 3 $\frac{1}{2}$ to 12 $\frac{1}{2}$

Insoles Ortocenter



Cod: 1845

Indications:

Relief in the pain of the heels.

Features:

Micro Foam and E.V.A.

Numbers from 3 ½ to 12 ½

Elastic Strap for Haluxx Valgo Traction



Cod: 1851

Indications:

Phalangian metatarsal bursitis of the large artelho, rheumatoid arthritis deformity, postoperative surgery of Haluxx Valgo.

Features:

Elastic, Micro Foam .

Sizes:

Small, Medium and Large

Gutter



Immobilization of the ankle and feet, replaces the short orthosis for positioning.

Features:

Foam coated polypropylene, Velcro closure.

Sizes:

Tailored

Short Immobilizer Boot Height 33 cm

Indications:



Cod: 555

Foam with duraluminium fins and rigid plastic, Velcro closure.

Ankle sprains, prevention for treatment of rheumatoid arthritis

and postoperative deformities.

Sizes:

Features:

Extra Small, Small, Medium, Large and Extra Large

Long Immobilizer Boot Height 39 cm



Indications:

Ankle sprains, prevention for treatment of rheumatoid arthritis and postoperative deformities.

Features:

Foam with duraluminium fins and rigid plastic, Velcro closure.

Extra Small, Small, Medium, Large and Extra Large

Long Immobilizer Boot Robocenter

Indications:



Cod: 558

Plastic injected Base, Foam with duraluminium fins and rigid plastic, Velcro closure.

Ankle sprains, prevention for

treatment of rheumatoid arthritis

and postoperative deformities.

Sizes:

Features:

Small, Medium, Large and Extra Large

Sling for Double Inguinal Hernia



Cod: 600

Indications:

Inguinal hernia support, avoiding growth.

Features:

Resistant elastic, polyurethane pellets lined soft and anatomymyas. Velcro closure.

Sizes:

80, 85, 90, 95, 100, 105 and 110 cm Circumference

Sling for Inquinal Hernia



Cod: 610

Indications:

Inguinal hernia support, avoiding growth.

Features:

Resistant elastic, polyurethane pellets lined soft and anatomymyas. Velcro closure.

Sizes:

80, 85, 90, 95, 100, 105 and 110 cm Circumference (Right and left)

Special Sling for Inguinal Hernia



Cod: 615

Indications:

Inguinal hernia support, avoiding growth.

Features:

Resistant elastic, Steel sheet, polyurethane pellets lined soft and anatomymyas. Velcro closure.

Sizes:

80, 85, 90, 95, 100, 105 and 110 cm Circumference (Right and left)

Sling for Umbilical Hernia



Cod: 620

Indications:

Hernia support in the navel region, avoiding growth.

Features:

Resistant elastic, polyurethane pellets lined soft and anatomymyas. Velcro closure.

Sizes:

80, 85, 90, 95, 100, 105 and 110 cm Circumference

Underwear for Inquinal Hernia



Cod: 625

Indications:

Inguinal hernia support, avoiding growth.

Features:

Coton, elastic, polyurethane pellets lined soft and anatomymyas.

Sizes:

Small, Medium, Large and Extra Large. (Right, left and double)

Scrotal Suspenders



Cod: 628

Indications:

Post-operative of scrotal hernia and vasectomy.

Features:

Cotton and adjustable elastic.

Sizes:

Small, Medium, Large and Extra Large.

Scrotal Suspenders - Without Regulating



Cod: 629

Indications:

Post-operative of scrotal hernia and vasectomy.

Features:

Cotton.

Sizes:

Small, Medium, Large and Extra Large.

ORTOCENTER

 $\hbox{\tt ORTHOPEDIC PRODUCTS}$

Ortocenter Aparelhos Ortopédicos Ltda. R. Baronesa de Bela Vista, 55 - São Paulo - SP - 04612-000 vendas@ortocenternet.com.br

Call: +55 (11) **5095 8555** www.ortocenternet.com.br

Knee Pad in Short Neoprene Smooth



Indications:

Instabilities and mild traumatic affections of the knee. treatment of arthropathies of the hemophiliac, prevention of recurrence in sports practice.

Features:

Neoprene Coated

Small, Medium, Large and Extra Large

Short Neoprene Knee Pad with Holes and Patellar Support



Instabilities of the Kneecap syndrome of Osgood Schlatter, treatment of Tendonitis, prevention of recurrence in sports practice.

Features:

Neoprene Coated, Compression pad and U-shaped support.

Small, Medium, Large and Extra Large

Long Neoprene Knee Pad with Bracket for the Kneecap



Indications:

Instabilities of the Kneecap, tendonitis, treatment of arthropathies of the hemophiliac, prevention of recurrence in sports practice.

Features:

Neoprene Coated, compression cushion and support for the kneecap.

Small, Medium, Large and Extra Large

Neoprene Hinged Knee Pad with Velcro Straps



Indications:

Traumatic knee disorders, chronic rotatory instabilities, arthrosis, arthritis and prevention of recurrence in sports practice.

Features:

Neoprene Coated, orifice, hinged made in PVC, Velcro closure.

Sizes:

Small, Medium, Large and Extra Large

Short Neoprene Knee Pad with Holes



Indications:

Instabilities of the Kneecap and mild traumatic affections of the knee, treatment of arthropathies of the hemophiliac, prevention of recurrence in sports practice.

Features:

Neoprene Coated

Small, Medium, Large and Extra Large

Short Neoprene Knee Pad with Bracket for the Kneecap



Indications:

Instabilities of the Kneecap, treatment of Tendonitis, prevention of recurrence in sports practice.

Features:

Neoprene Coated, Compression pad and support in E.V.A. for the kneecap

Small, Medium, Large and Extra Large

Long Neoprene Knee Pad with Holes



Indications:

Instabilities of the Kneecap and mild traumatic affections of the knee, treatment of arthropathies of the hemophiliac, prevention of recurrence in sports practice.

Features:

Neoprene Coated.

Small, Medium, Large and Extra Large

Long Neoprene Knee Pad with Holes and Bracket for the Kneecap



Indications:

Instabilities of the Kneecap, treatment of Tendonitis, prevention of recurrence in sports practice.

Features:

Neoprene Coated, Compression pad and support in E.V.A. for the kneecap.

Small, Medium, Large and Extra Large

Neoprene Coxal



Cod: 1650

Indications:

Bruises, muscle distention. Warms the musculature, prevention of recurrence in sports practice.

Features:

Neoprene Coated.

Sizes:

Small, Medium, Large and Extra Large

Long Neoprene Knee



Indications

Instabilities of the Kneecap and mild traumatic affections of the knee, treatment of arthropathies of the hemophiliac, prevention of recurrence in sports practice.

Features:

Neoprene Coated.

Sizes:

Small, Medium, Large and Extra Large

Long Neoprene Anklet



Cod: 1685

Indications:

Compression of the ankle joint in the treatment and prevention of edema of sprains, tendinitis, bursitis and rheumatoid arthritis, used in sports activities and prevention of recurrence.

Features:

Made of Neoprene Coated.

Sizes:

Small, Medium, Large and Extra Large

Short Neoprene Anklet



Cod: 1680

Indications:

Compression of the ankle joint in the treatment and prevention of edema of sprains, tendinitis, bursitis and rheumatoid arthritis, used in sports activities and prevention of recurrence.

Features:

Made of Neoprene Coated.

Sizes:

Small, Medium, Large and Extra Large

Long Neoprene Knee Pad with Holes and U Pattelar Suport



Cod: 1691

Indications:

Instabilities of the Kneecap, treatment of Tendonitis, prevention of recurrence in sports practice.

Features:

Neoprene Coated, Compression pad and support in E.V.A. for the kneecap.

Sizes

Small, Medium, Large and Extra Large

Neoprene Hinged Knee Pad with Cross Straps



Cod: 1695

Indications

Traumatic knee disorders, chronic rotatory instabilities, arthrosis, arthritis Kneecap stabilizing, lateral and medial fin in plastic and prevention of recurrence in sports practice.

Features:

Neoprene Coated, orifice, cross straps and Velcro closure.

Sizes:

Small, Medium, Large and Extra Large

Neoprene Coxal



Cod: 1698

Indications:

contusions, bursitis, arthritis and other affections of the elbow, prevention of recurrence in sports practice.

Features:

Neoprene Coated.

Sizes

Small, Medium, Large and Extra Large

Neoprene Wrist Protector With Thumb Guard



Cod: 1700

Indications:

Mild cuff affections, prevention of recurrence in sports practice.

Features:

Neoprene Coated, Velcro closure.

Sizes:

Single Size (Right and left)

Neoprene Wrist Protector with PVC Backing on Thumb



Indications:

Mild cuff affections, prevention of recurrence in sports practice.

Features:

Neoprene Coated, PVC backing on Thumb Velcro closure.

Sizes:

Single Size (Right and left)

Neoprene Wrist Without Thumb



Indications:

Mild cuff affections, prevention of recurrence in sports practice.

Features:

Neoprene Coated, Velcro closure.

Sizes:

Single Size (Right and left)

Simple Neoprene Wrist



Indications:

Mild cuff affections, prevention of recurrence in sports practice.

Features:

Neoprene Coated, Velcro closure.

Sizes:

Single Size

Strap for Tennis Elbow in Neoprene



Indications:

Epicondylites and other affections the elbow, prevention of recurrence of sports practice.

Features:

Neoprene coated, Velcro closure. Bilateral.

Sizes:

Small and Large

Immobilizer Splint in Neoprene for Handle



Cod: 1715

Indications:

Tenosynovitis for wrist, arthroses, ligamentous lesions, corrections of the vicious positions in rheumatoid arthritis, region affections and prevention of recurrence in the return of activities.

Made of Tissue, Velcro closure and internal reinforcement in rigid plastic.

Sizes:

Small, Medium, Large and Extra large. (Right and left)

Neoprene Glove with Thumb



Cod: 1735

Indications:

Prevention of recurrence in sports practice.

Features:

Neoprene Coated, With leather stand on the palm of the hand and thumb, Velcro closure.

Sizes:

Small, Medium, Large (Right and left)

Neoprene Knee Pad with Tensioner



Indications:

Traumatic knee disorders, chronic rotatory instabilities, arthrosis, arthritis and prevention of recurrence in sports practice.

Features:

Neoprene Coated, with tensioner.

Small, Medium, Large and Extra Large

Neoprene Lumbar Range



Indications: Lumbar spine Stabilization, post

partum and post surgical.

Features:

Made of Neoprene coated, velcro closure.

Sizes:

Small, Medium, Large and Extra large

Cervical Collar Schantz or Thomas



Cod: 010

Indications:

Traumatisms, torticollis, arthritis, arthroses and other cervical spine affections.

Features:

Made of Sturdy flexible plastic, upholstered on the edges, support cushion on the external, height adjustment and velcro closure.

Sizes:

Extra Small

Cervical Collar with Mentonian Support



Cod: 015

Indications:

Traumatisms, torticollis, arthritis, including rheumatoid, arthroses and other cervical spine affections

Features:

Made of Sturdy flexible plastic, upholstered on the edges, Mentonian support to reduce movement, height adjustment and velcro closure.

Sizes:

Extra Small

Structural Necklace



Cod: 040

Indications:

Support for the transport situation, tours and learning activities.

Features:

Made of Foam lined with tissue and velcro closure.

Sizes:

Single Size (Child)

Sling for Umbilical Hernia



Cod: 085

Indications:

Hernia child support in the navel region, avoiding growth.

Features:

Resistant elastic, pellet lined soft and anatomymyas. Velcro closure.

Sizes:

Small, Medium and Large (Child)

Velpeau type Immobilizer in Upholstered Nylon (Infant)



Cod: 201

Indications:

Immobilization of the upper limb, contusions dislocation, fractures of the humerus, elbow, bursitis and sequelae of the A.V.C.

Features:

Made of nylon with foam, plastic buckle with one-handed manipulation, Velcro closure, bilateral.

Sizes:

Single Size – Child

Frejka Diaper



Cod: 400

Indications:

Disluxations and congenital subluxations of the hip, preventive use in the first months of life.

Foatures:

Cotton Tissue, and foam inner cushion. Velcro closure.

Sizes:

Small, Medium, Large and Extra large

Pavlik Suspenpensory



Cod: 402

Indications:

Disluxations and congenital subluxations of the hip, contributes to the conformation of the acetabulum and reduces the osteochondritis risk.

Features:

Cotton Tissue. Velcro closure.

Sizes:

Small, Medium, Large and Extra large

Elastic Sling Handset for Lower Limbs Rotation



Cod: 405

Indications:

Correction of the internal rotation of the lower limbs through reeducation of posture. Daytime and nighttime use.

Features:

Resistant Elastic, ajustable. Velcro closure.

Sizes:

Small, Medium, Large and Extra large

∞

Emergency Cervical Collar - Rescue



Cod: 037

Indications:

Support for the cervical spine during r emoval of the injured traumatisms.

Features:

Made of High density polyethylene, coated in E.V.A. and Mentonian support. With wide frontal opening for analysis of the choroid pulse, access to trachea and velcro closure.

Sizes:

Extra Small, Small, Medium, Large and Extra Large

Long Plank for Emergency



Cod: 2023

Indications:

Transport of accidents during rescue.


Features:

Made of Naval plywood
Dimensions 1.85 x 0, 46m,
rounded outer edges. Waterproofed
with holes for support and passage
of the ground contact support straps
and passage for straps.
Accompanies 3 fastening belts.

Sizes:

1.85 x 0. 46m

Plank for Emergency - Child



Cod: 2024-03

Indications:

Transport of accidents during rescue.

Features:

Made of Polyethylene Dimensions 1.28 x 0, 44m, r ounded outer edges. Holes for support and passage of the ground contact support straps and passage for straps. Accompanies 3 fastening belts.

Sizes:

1.85 x 0, 46m

Belt for Fastening on Plank Removal



Cod: 2034

Indications:

Against accidental loosening and to ensure prompt and immediate loosening by the paramedics.

Features:

Set of 3 belts made of highly resistant tapes with coupling system (male and female) that supports up to 50 kg, easy to regulate and translucent to RX.

Dorsal Immobilizer Vest - KED



Cod: 2030

Indications

Indicated in the immobilization of the spine in the rescue and rapid, safe and efficient extraction of the injured.

Features:

Sturdy nylon, and stiffened with natural wood allowing for greater product durability. Headrest, translucent product to RX, Accompanies tote bag with zippers.

Sizes:

Adult

Long Plank for Emergency Polyethylene



Cod: 2025

Indications:

Transport of accidents during rescue.

Features:

Made of Polyethylene, Dimensions 1.88 x 0, 44m, rounded outer edges. Waterproofed with holes for support and passage of the ground contact support straps and passage for straps. Accompanies 3 fastening belts.

Sizes:

1.88 x 0, 44m

Dorsal Immobilizer Vest - KED Child



Cod: 2032

Indications:

Indicated in the immobilization of the spine in the rescue and rapid, safe and efficient extraction of the injured.

Features:

Sturdy nylon, and stiffened with natural wood allowing for greater product durability. Headrest, translucent product to RX, Accompanies tote bag with zippers.

Siz

Sizes: Child

Emergency Farmers Market



Cod: 2040

Indications:

Indicated for major events in places of great agglomeration of people.

eatures:

Resistant nylon with impact insulators, has 2 divisions. Used for first aid.

Sizes:

90 x 28 x 37 cm

Emergency Bag - Slin



Cod: 2040-1

Indications:

Indicated for major events in places of great agglomeration

Features:

Resistant nylon with impact insulators, has 2 divisions. Used for first aid.

Sizes:

90 x 25 x 35 cm

Emergency Open Bag



Cod: 2040-A

Indications:

Indicated for major events in places of great agglomeration of people.

Features:

Resistant nylon with impact i nsulators, has 2 divisions. Used for first aid.

Emergency Bag - Slin



Cod: 2040-B

Indications:

Indicated for major events in places of great agglomeration of people.

Resistant nylon with impact insulators, has 2 divisions. Used for first aid.

Backpack Pre-Hospital Care



Cod: 2041

Indications:

Used in the pre-hospital Emergency care and by removal companies.

Features:

Resistant nylon with impact insulators, has PVC Pockets.

MOCHILA de Medicamento



Indicações:

Utilizado para o transporte de ampolagens e medicamento em geral.

Características:

Confeccionado em nylon resistente, revestida com isolante à impactos e possuem bolsos em P.V.C.

Cod: 2042

Trauma Tote Bag



Indications:

Features: Resistant nylon with impact

insulators, has 1 division.

immobilization KIT and procedure.

Suitable for transport of

Cod: 2050

Trauma Tote Open Bag



Indications:

Suitable for transport of immobilization KIT and procedure.

Resistant nylon with impact insulators, with several sub divisions.

Cod: 2050-A

Emergency Bag - 192



Indications:

Used in the pre-hospital Emergency care and by removal companies.

Features:

Resistant nylon with impact insulators, with five internal partitions and supporting feet. Strap with 40 cm.

Cod: 2041-2

Silicone Finger Separators



Cód.: **1837**

Indications:

Separate from crooked and curled fingers and fix bunions.

Features:

Made of 100% silicone.

Sizes:

Small, Medium, Large

Silicone Heel without Edge



Cod: 1840

Indications:

Relief of Calcaneal Spur, bursitis and plantar fascitis.

Features:

Made of 100% silicone, does not change the thickness and does not heat.

Sizes:

Small, Medium, Large

Silicone Insole with Soft Bow



Cod: 1841

Indications:

Impact absorption, better comfort for knee and spine pain.

Features:

Made of 100% silicone, does not change the thickness and does not heat.

Sizes:

Numbers from 3 ½ to 12 ½.

Silicone Flat Insole



Cod: 1842

Indications:

Impact absorption, better comfort for knee and spine pain. Suitable for sports practice.

Features:

Made of 100% silicone, does not change the thickness and does not heat.

Sizes:

Numbers from 3 $\frac{1}{2}$ to 12 $\frac{1}{2}$

Silicone Heel



Cod: 1846

Indications:

Impact absorption, Bursitis, Achilles tendinitis, relief of the Calcaneal Spur. Suitable for sports practice.

Features:

Made of 100% silicone, does not change the thickness and does not heat.

Sizes:

Small, Medium, Large and Extra Large

Silicone Heel



Cod: 1847

Indications:

Impact absorption, Bursitis, Achilles tendinitis, relief of the Calcaneal Spur. Suitable for sports practice. with blue Dot.

Features:

Made of 100% silicone, does not change the thickness and does not heat.

Sizes:

Small, Medium, Large and Extra Large



Ortocenter Aparelhos Ortopédicos Ltda.

Rua Baronesa de Bela Vista, 55 - São Paulo - SP CEP 04612-000 - Indústria Brasileira

Call: +55 11 **5095 8555 - www.ortocenternet.com.br** vendas@ortocenternet.com.br



NOVELTY NEOPRENE



Abdominal Brace for Obese and Post Surgical



Cod: 066

Indications:

Postoperative, morbid obesity, stomach reduction surgery.

Features:

Made of elastic with Resistant elastics, smooth abdominal compression, easy on-site inspection and hygiene, Velcro closure.

Sizes:

Small, Medium, Large and Extra large

Lumbar Brace with Suspenders for Safety at Work in Elastic



Cod: 067 / 068

Indications:

Spine protection prevents injuries caused by excessive physical exertion.

Features:

Made of elastic with Resistant elastics, Fasteners and velcro adjustment and adjustable suspenders.

Sizes

Small, Medium, Large and Extra large

Lumbar Brace with Suspenders for Safety at Work in Canvas



Cod: 069

Indications:

Spine protection prevents injuries caused by excessive physical exertion.

Features:

Made of Heavy Duty Canvas fasteners and velcro adjustment and adjustable suspenders.

Sizes:

Small, Medium, Large and Extra large.

Brace for Neoprene Therapy



Cod: 1745

Indications:

External use in the treatment of muscle spasms luxations, recent contusions and physiotherapy.

Features:

Neoprene coated with pouch for hot or cold water. Velcro closure.

Sizes:

Single Size.

Bath Collar for Tracheostomy Surgery



Cod: 042

Indications:

During bathing for patients who underwent tracheostomy surgery.

Features:

PVC 2.5 mm with polyester lining a nd velcro closure,

Sizes:

Small and Large.

Anti Snoring Track



Cod: 076

Indications:

During sleep for those who snore.

Features:

Neoprene Coated.

Sizes:

Single Size

Immobilizer Cuff and Thumb Splint Neoprene



Cod: 1717

Indications:

Tenosynovitis for wrist, arthroses, ligamentous lesions, rheumatoid arthritis, region affections.

Features:

Made of Neoprene coated, Velcro closure and internal reinforcement in rigid plastic.

Sizes:

Small, Medium, Large and Extra large. (Right and left)

Strap for Pregnant Woman with Lumbar Reinforcement



Cod: 116

Indications:

Abdomen support, relief of low back pain, relief of abdominal muscles and sp+. Posture Improvement.

Features

Made with Front support in solid tissue with elastic rods,

Sizes:

Single Size

Elbow Protector



Cod: 741

Indications:

Post surgical and avoids escaras, pressure ulcer.

Features:

Made with post surgical and avoids escaras, pressure ulcer. Velcro Closure.

Sizes:

Single Size

Heel Protector



Cod: 743

Indications:

Post surgical and avoids escaras, pressure ulcer.


Features:

Made with post surgical and avoids escaras, pressure ulcer.
Velcro Closure.

Sizes:

Single Size

Mattress Type Egg box



Cod: 5070

Indications:

Avoids escaras, pressure ulcer.

Features:

Soft Foam.

Sizes:

Single Size

Hip Abduction Pad



Cod: 745

Indications:

Cushion for positioning the lower limbs after surgical hip arthroplasty, prevention of skin ulcerations caused by pressure.

Features:

Made with Polyurethane foam and adherent fasteners.

Sizes:

Height – 170mm Width – 220mm Length – 310mm

Pro Splint Type Immobilizer Splice



Cod: 2016

Indications:

For adults and children for use in knee, leg, ankle, arm, wrist, and hand. Easy application, lightweight, washable, Velcro closure.

Features:

Pró Splint.

Sizes:

6-Piece KIT and 10-Piece KIT

Bracelet and Anklet for Restriction of Bedridden Patient



Cod: 317

Indications:

Fix patient in bed.

Features:

Canvas in bracelet shape and velcro closure.

Sizes:

Adult and Child

Disposable Plus Restrictor



Cod: 319-03

Indications:

Fix patient in bed.

Features:

5mm polyurethane foam and linen tissue with cotton strips.

Sizes

Adult 35 x 10 x 200 cm; Child 25 x 7 x 150 cm; Newborn 18 x 4,5 x 100 cm

ORTOCENTER

ORTHOPEDIC PRODUCTS

Ortocenter Aparelhos Ortopédicos Ltda.

R. Baronesa de Bela Vista, 55 - São Paulo - SP - 04612-000 vendas@ortocenternet.com.br

Call: +55 11 5095 8555 www.ortocenternet.com.br

Immobilizer Head



Cod: 2010

Indications:

For immobilization of the head, cervical spine and transport on a rescue board.

Features:

Made of Soft Foam, Velcro closure.

Sizes:

Adult, Child

Splint for Provisional Immobilization with E.V.A.



Cod: 3030

Indications:

Splint for provisional immobilization.

Features:

Made of wire and coated with E.V.A.

Sizes:

Extra Small – 30 x 8cm; Small – 43 x 8cm; Medium – 63 x 9cm; Large – 86,5 x 10cm Extra Large 102 x 11,5; Extra Extra Large – 120 x 12cm

Foam Splint for Provisional Immobilization



Cod: 3050

Indications:

Splint for provisional immobilization.

Features:

Made of foam.

Sizes:

Extra Small – 30 x 8cm; Small – 43 x 8cm; Medium – 63 x 9cm; Large – 86,5 x 10cm Extra Large 102 x 11,5; Extra Extra Large – 120 x 12cm

Cover and Immobilization KIT



Cod: 2057

Indications:

Cover and immobilization KIT for guarding (immobilizer head plank and rescue collars).

Features:

Made of Sturdy canvas with zipper closure.

Sizes:

Single size

KIT Childbirth



Cod: 100105

Indications:

Case with absorbent, apron, scalpel, pair of surgical gloves, ten sterile compresses, umbilical fastener, alcohol sachet and purge bag.

All materials are disposable.

Features:

Materials are disposable.

Aluminized Thermal Blanket



Cod: 2015

Indications:

Keep the body warm.

Features:

Metallized and aluminium-coated bioriented polypropylene film.

Sizes:

2,10 x 1,40m 25 Microns Thickness

Triangular Bandage



Cod: 2017

Indications:

Provisional immobilization.

Features:

Cotton tissue.

Sizes:

2,00 x 1,40 x 1,40m; 1,40 x 1,00 x 1,00m; 1,00 x 0,70 x 0,70m

KIT Rescuer



Cod: 2038

Indications:

Fanny pack adapted with belts.

Features:

1 Pair of gloves, 1 mask.

Size:

Single Size

Belts for Spider type Rescue Board



Indications:

For immobilization of the body in rescue board.

Features:

Made of Highly resistant Tapes. Velcro closure.

Sizes:

Adult, Child

Belts for Long Rescue Board



Cod: 2035

Indications:

For immobilization of the body in long rescue board.

Features:

Made of Highly resistant Tapes. Velcro closure.

Sizes:

Adult

THE ORTOCENTER SAVING LIVES





The ORTOCENTER offers:

- · Customer Engagement.
- · Customer Engagement.
- · Telesales for Customer.
- · Economy in the purchase of products.
- · Products with factory technology.
- · Delivery with deadlines and speed.
- · Agility in products on request.
- · Payment Deadlines.
- · Customer Technical Support.
- · Internet and social networks to streamline contacts.

Call: +55 11 **5095-8555**

The Ortocenter caring for You

The Ortocenter caring for You



ORTOCENTER ORTHOPEDIC PRODUCTS

Call: +55 (11) 5095 8555

www.ortocenternet.com.br - Parceiro Virtual: www.ortoplusnet.com.br vendas@ortocenternet.com.br - vendas@ortoplusnet.com.br























Ortocenter Aparelhos Ortopédicos Ltda.

Rua Baronesa de Bela Vista, 55 - São Paulo - SP CEP. 04612-000 - Indústria Brasileira

